

Kinetics Proposal 1

Your name: Cieran Wong

Your email: cieranwong@arizona.edu

Your lab partner's name(s): Nick Samuels
nssamuels014@arizona.edu

Your lab partner's email(s):

Your lab instructor's name: Pefro

Your lab section: 2A

*All work must be **very neat and organized**. If you need to collect your thoughts, please use a separate sheet of paper. Proposals are a **group effort**. Please submit the completed document as a PDF to the **Kinetics Proposal 1 D2L DropBox** folder before the scheduled end of lab.*

1. In a complete, well-written sentence, summarize in your own words the **overall goal(s)** for the *Kinetics Project*.

Overall goal is to determine the shelf life of two drugs in an alkaline based solution at room temperature by finding out the rate order and observed rate constant. This information will help us decide how to best preserve the drug so that it does not degrade through changing the concentration of the alkaline solution and the concentration of the drugs used.

Do drugs

2. In your own words, the **goal for this first session** of the *Kinetics Project* is...

Goals for the first session are to find the NaOH concentration and volume ratio of Phenolax/NaOH and Pyoctanin/NaOH that give initial absorbance values between 1-1.1 and degrades over a period of 150-250 seconds to an absorbance value of 0.2-0.1. Then from the results, determine a proposal to find the rate order, observed rate constant, and shelf life.

3. **Semi-Quantitative Exploration Results.** Clearly summarize the **key results** from your systematic semi-quantitative exploration of 1) at least 3 different *volume ratios of the Phenolax and NaOH solutions/dilutions* on the Phenolax degradation rate, and 2) at least 3 different *volume ratios of the Pyoctanin and NaOH solutions/dilutions* on the Pyoctanin degradation rate. Be sure to specify the diluted NaOH concentrations indicating the volumes used to prepare the NaOH dilutions - there should be at least two (2) different NaOH dilutions explored for Pyoctanin and Phenolax. Hint: What happens with the observed color intensity over time for different ratios and NaOH dilutions? What volume of Pyoctanin and Phenolax in relation to the NaOH concentration and volume are within the appropriate absorbance range to collect your data?

Table 1: Phenolax and Pyoctanin NaOH Degradation Semi-Quantitative Exploration Data and Observations.

NaOH Concentration	Volume of NaOH	Volume of NanoPure	Volume of Phenolax	Volume of Pyoctanin	Color at time = 0 seconds	Color at time = 250 seconds
1.0M	3.00 mL	0.00 mL	29	64	Pheno: Bright pink magenta color Pyoc: Dark Indigo purple blue	Color is gone clear once like it was before the adding of the drug

						Color was gone at around 180-190 seconds
0.75M	2.25 mL	0.75 mL	29	64	Pheno: Bright pink Pyoc: Dark Indigo purple blue	Color fully cleared out around 190 seconds
0.5M	1.50 mL	1.50 ml	29	64	Pheno: Bright pink Pyoc: Dark Indigo purple blue	For Pheno at 0.5 M there was still color after the 250 seconds we observed for. Pyoc all colored cleared out before the 250 seconds ended fully gone around 180 seconds. The 0.5 molarity ratio took the longest for color to disappear especially when we looked at color change for pheno

To calculate the volume of Phenolax needed:

$$\text{Phenolax} = \frac{0.75}{5.02 \times 10^4 \text{ cm}^{-1} \text{ M}^{-1} \times 1.00 \text{ cm}} = 1.49 \times 10^{-5} \text{ M} = \frac{1.49 \times 10^{-5} \text{ M} \times 3 \text{ mL}}{1.57 \times 10^{-3} \text{ M}} = 0.02847 \text{ mL} \times 1000 = 28.47 \mu\text{L}$$

To calculate the volume of Pyocyanin:

$$\text{Pycotanin} = \frac{0.75}{8.65 \times 10^5 \text{ cm}^{-1} \text{ M}^{-1} \times 1.00 \text{ cm}} = 8.6705 \times 10^{-7} \text{ M} = \frac{8.6705 \times 10^{-7} \text{ M} \times 3 \text{ mL}}{4.08 \times 10^{-5} \text{ M}} = 0.06375 \text{ mL} \times 1000$$
$$= 63.75 \mu\text{L}$$

4. **Proposal 1.** Based on your semi-quantitative exploration, **propose a plan**, and **justify each step** to experimentally determine n , k_{obs} and t_{90} at room temperature for **PHENOLAX**, where the condition $[NaOH] \gg [Phenolax]$ exists as described in the *Kinetics Guide*. **Please NUMBER your procedural steps.**

Procedural Step	Justification based on data/observations, or technical instructions, or conceptual understanding
<ol style="list-style-type: none"> 1) Fill a 20 ml vial with NaOH 2) Gather 3, 10 ml vials. 3) Choose 3 Molarity ratios, e.x. 1.0 M, 0.75 M, 0.5 M. 4) Calculate the appropriate volumes of the systems for a total volume of 3ml using $C_1 \cdot V_1 = C_2 \cdot V_2$ 5) For Phenolax, set aside three vials and fill it with the three different molarity ratios. 6) Use Beer's Law to calculate the molarity of Phenolax solutions based on the chosen absorbance levels 7) Use the molarity of the Phenolax to calculate the volume of the respective solutions using the following formula: $C_d V_d = C_s V_s$ 8) Take the molarity in step 5, times that by the total volume of the NaOH and Nanopure solution and divide it by the molar mass of the respective solutions. 9) In the three vials from step3 containing the NaOH and Nanopure mixture, add in the determined volume of Phenolax in Step 8. 10) Add the solution to a cuvette and place it in the spectrometer, recording data for absorbance vs. time. Set a stopwatch for 250 seconds, observing any changes in color. 11) Using the observed absorbance value as it degrades over time, we can then calculate the rate order, the K_c values, and the shelf life. 	<ol style="list-style-type: none"> 1) Gathering an appropriate amount of NaOH for the upcoming experiment 2) These vials are used to create our 3ml systems in which we will place our drug (phenolax) in 3 different molarity ratios 3) In order to examine how different ratios of phenolax to NaOH affect different the degradation of the drug 4) You are to determine the exact volumes for each corresponding ratio so that your ratios have a small margin of error, and readings will be accurate 5) Clearly separating the ratio systems from each other ensuring not to have them mixed up 6) Beer's Law gives us the exact molarity of our pheno solution at a specific absorbance level 7) Using the formula of: $C_d V_d = C_s V_s$ you can determine the exact volume of Phenolax. 8) This will give us the volume in milliliters, multiply by 1000 to get the volume in microliters. 9) This will allow us to observe color changes as it will provide a visualization of Phenolax degrading over time. 10) This will allow us to record the absorbance value from the instant Phenolax is placed into the NaOH solution. It needs to happen instantly as Phenolax can degrade fairly quickly. 250 seconds should be ample time to fully see the Phenolax in the NaOH solution degrade. 11) The rate order, K_c values and the shelf life are the values we are looking for and the absorbance level from the spectrometer can help us determine it. We will use the following equations: $t_{90} = \frac{[X]_0}{10k_{obs}}$ $t_{90} = -\frac{\ln(0.9)}{k_{obs}}$

$$t_{90} = \frac{0.111}{k_{\text{obs}}[X]_0}$$

The three equations will give us three different graphs that we can use to find the R² value. This value will help us determine which formula is applicable as we can see which R² value is closest to 1.